



**electric driving support**  
on all patient lifts



become exclusive    
**client** or **partner**

# intuitive **touch&move** driving support

**touch&move** is a unique, intuitive electrically powered drive that makes working with a patient lift easier and more pleasant. The push bar of a patient lift with **touch&move** driving support has sensors that measure the magnitude and direction of the force exerted on it. A gentle push against the bar is enough to set the additional wheel under the patient lift in motion, either forwards, backwards or sideways. That may sound technical but using it is very intuitive.

Watch our models  
and versions here



## THE BENEFITS OF TOUCH&MOVE

- Makes manoeuvring lighter and easier
- Reduces physical strain injuries
- Increases efficiency of caregivers
- Intuitive through (years of) continued development
- Contributes to healthy working

## **S**MOVE

The **S**MOVE is an active patient lift (standing lift), providing the patient with comfortable support to get them from a seated to an upright position. Encouraging their own capabilities and the standing function are important aspects. The **S**MOVE has electrically adjustable legs that can be spread apart and a high/low function for the lifting arm.

MAXIMUM LIFTING WEIGHT **200 kg**



## **Gentle**MOVE

The **Gentle**MOVE is a passive patient lift for moving patients safely and comfortably. **Gentle**MOVE's lifting arm travels in an arc. This rotating movement gives you an even greater range when the patient lift is driven from the front. This allows the patient to be positioned right at the back of the chair. The spread of the **Gentle**MOVE's legs is electrically adjustable, as are the high/low functions of the lifting arm.

MAXIMUM LIFTING WEIGHT **200 kg**



## **Power**MOVE

The **Power**MOVE is a passive patient lift for moving patients safely and comfortably. The mast allows a translational movement of the lifting arm. The patient is always lifted the same distance from the mast, regardless of height, and they can therefore be positioned very precisely. The **Power**MOVE requires little maintenance and the spread of its legs is electrically adjustable, as are the high/low and tilting functions of the arm and spreader bar.

MAXIMUM LIFTING WEIGHT **200 kg**



## **Power**MOVE *Mini*

The **Power**MOVE *Mini* is a passive column patient lift. The **Power**MOVE *Mini* is more compact than the **Power**MOVE, thanks to its shorter arm and legs. This makes the **Power**MOVE *Mini* highly manoeuvrable – ideal for smaller spaces. The mast allows a translational movement of the lifting arm. The patient is always lifted the same distance from the mast, regardless of height, and they can therefore be positioned very precisely.

MAXIMUM LIFTING WEIGHT **150 kg**





# We are JOYinCARE

**JOYinCARE** is all about making care more enjoyable. We do this by tackling physical strain step by step, from the source. **Why?** To keep care sustainable and affordable, both now and in the future, while maintaining as much self-esteem as possible for the patient. Together with caregivers, patients and paramedics, we devise solutions for each unique need. We perform our work with a mission in mind:

JOYinCARE, for sustainable care.

become exclusive    
**client** or **partner**

[www.joyincare.com](http://www.joyincare.com)

Industrieweg 47, 8071 CS Nunspeet, The Netherlands

✉ [info@joyincare.com](mailto:info@joyincare.com)

☎ +31 (0)341 - 76 80 80

 Made in Holland